














Overview: Supported inkjet printers (1)

Supported printers:	MI 9330	MI 9029	MI 9410	MI 9450
<ul style="list-style-type: none"> MI 9029 MI 9330 MI 9410 MI 9450 MI 9018 (discont.) MI 9020 (discont.) MI 9028 (discont.) MI 9030 (discont.) MI 9232 (discont.) MI 9341 (discont.) 	 www.markem-ijaje.com/	 www.markem-ijaje.com/	 www.markem-ijaje.com/	 www.markem-ijaje.com/
Performance	4 Lines Alphanumeric wide choice of characters Graphics/Logo printing Side wall printing	4 Lines Alphanumeric wide choice of characters Graphics/Logo printing Side wall printing	5 Lines Alphanumeric, bar codes, 2D codes Graphics/Logo printing Side wall printing	5 Lines Alphanumeric, bar codes, 2D codes Graphics/Logo printing Side wall printing
Other features	Bench top Touch screen USB port RS-232 / Ethernet	Bench top Key pad USB port RS-232	Bench top Touch screen USB and SD port Ethernet	Bench top Touch screen USB and SD port Ethernet
Dimensions (HxWxD)	595 x 303 x 280 mm	595 x 307 x ??? mm	555 x 360 x 357 mm	555 x 360 x 357 mm
Weight (kg)	18.5	18.5	25	25
List Price (EUR)	<10'000	<10'000	>10'000	>10'000
Maintenance	ink exchange in < 6 min ink runs > 8000 hours	ink exchange in < 6 min ink runs > 8000 hours automatic cleaning system	Maintenance: every 18 month Head cleaning recommended every second week Head can be left attached on MEDIAJET	Maintenance: every 18 month Head cleaning recommended every second week Head can be left attached on MEDIAJET
INTEGRA comments	Recommended MI printer, lower consumable and maintenance costs. Industrial interface must be ordered.	Industrial interface must be ordered.	Industrial interface must be ordered.	Industrial interface must be ordered. Higher line speed and storage capacity.

Overview: Supported inkjet printers (2)

Supported printers:	Linx 89xx (89xx: -00, -10, -20, -40)	Linx 10	Domino Ax150i	Domino Ax350i
<ul style="list-style-type: none"> Linx 4900, Linx 5900 (printhead MK7) Linx CJ400, Linx 89xx (printhead MK11) Linx 10 (printhead MK13) Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i 	 www.linxglobal.com	 www.linxglobal.com	 www.domino-printing.com	 www.domino-printing.com
Performance	6 Lines Alphanumeric, barcodes, 2D codes Graphics/Logo 25 px Side wall / bottom printing	3 lines Code height 1.8 – 8.8 mm Upper/lower case, graphics Side wall / bottom printing	3 Lines Alphanumeric, bar codes, 2D codes Graphics/Logo printing Side wall printing	5 Lines Alphanumeric, bar codes, 2D codes Graphics/Logo printing Side wall printing
Other features	Benchtop Touch screen USB port Ethernet	Benchtop Touch screen Remote connectivity option Quick and easy to set up	Bench top Touch screen Optional: USB port, Ethernet, RS232, GPIO	Bench top Touch screen Optional: USB port, Ethernet, RS232, GPIO
Dimensions (HxWxD)	503 x 440 x 355 mm	228 x 551 x 189 mm	373 x 400 x 393 mm	411 x 430 x 331 mm
Weight (kg)	24	11	16	22
List Price (EUR)	>10'000	< 8'000	8'300	>10'000
Maintenance	Maintenance: every 24 months Head cleaning every 3 months	Service module (Easi-Change) to be exchanged after 18 months or 4.000 h / self-service possible	after 3500 working hours no external service necessary, 10 min. maintenance Head can be left attached on MEDIAJET	after 3500 working hours no external service necessary, 10 min. maintenance Head can be left attached on MEDIAJET
INTEGRA comments	Recommended, if 2D barcodes are needed. Reliable.	Recommended Linx printer Portable, lightweight, small footprint.	Recommended Domino printer Recommended, if bar codes or 2D codes needed.	Recommended, if bar codes or 2D codes needed.

Ordering information MEDIAJET inkjet printer connections kits

Supported printers:	Markem-Imaje (MI) 9018-9030, 9330	Markem-Imaje (MI) 9232/94xx	Linx side wall printing (with print head MK7, MK11, MK13)	Linx bottom printing (with print head MK7, MK11, MK13)	Domino
<ul style="list-style-type: none"> MI 9018, MI 9020, MI 9028, MI 9029, MI 9030, MI9330 MI 9232, MI 9410, MI 9450 Linx 4900, Linx 5900, Linx 89xx (print head MK7) Linx CJ400 (print head MK11) Domino A320i, Domino A420i, Domino Ax150i, Domino Ax350i 					
Part number					
Printer connection kit: with printer head fixing device and interface cable 	103080	103091	113845	113840	113060
Operating instructions	Operating instructions for connection kits are included and also available in the download center on www.integra-biosciences.com				
	103983	103980	103956	103981	103982
Cart (W x D x H (cm): 70 x 50 x 65-90) 	-	-	-	113841	-
INTEGRA comments	Compatible with: Ø 90 mm, Ø 60 mm or Ø 35 mm Petri dishes.	Compatible with: Ø 90 mm or Ø 60 mm Petri dishes.	Compatible with: Ø 90 mm or Ø 60 mm Petri dishes.	Compatible with: Ø 90 mm, Ø 60 mm or Ø 35 mm Petri dishes.	Compatible with: Ø 90 mm or Ø 60 mm Petri dishes.
Compatibility with MEDIAJET (manufactured)	After Dec 2005 / Feb 2007; (SN above: MJ : 10305000 / MJ vario : 11307000)	After Feb 2013; (SN above: MJ : 10313568 / MJ vario : 11313247)	After Feb 2013; (SN above: MJ : 10313568 / MJ vario : 11313247)	After Mai 2010; (SN above: MJ : 10310371 / MJ vario : 11310128)	After Feb 2013; (SN above: MJ : 10313568 / MJ vario : 11313247)